



Let Our Image Enhance Your Image

Superior image quality is the best way to describe this new MedRx Video Otoscope System

This user friendly system consists of a well balanced camera and probe that require no focusing. The miniature color video camera is light weight and has a precision white balance built in. The deluxe probe consists of an aluminum body and durable stainless steel conical tip that tapers down to 3mm.

The light source utilizes a 21 volt, 150 watt halogen bulb with an approximate life of 200 hours.



MedRx INCTM

www.medrx-int.com

MedRx International Sales
Salzufer 13/14, 10587 Berlin, Germany
Tel.: +49 30 70 71 46-50, Fax: -99
E-mail: medrx-sales@maico.biz