



Diagnostic & Testing Instrumentation

COMPUTER REQUIREMENTS

MedRx **Minimum** Computer Specs:

- Windows®-PC compatible computer
- Intel™ Dual Core, 1.8 GHz or better
- 2 GB RAM
- 5 GB free hard drive space
- Available 2.0 USB ports (2)
- CD-ROM or DVD-ROM drive
- Windows 7 or 8 Professional (32 or 64-bit)

MedRx **Recommended** Computer Specs:

- Windows®-PC compatible computer
- Intel™ Dual Core, 2.0 GHz or better
- 3 GB RAM
- 20 GB or more free hard drive space
- Available 2.0 USB ports (4)
- Graphics Adapter with Dedicated Video Memory
- DVD-ROM Drive
- High Speed Internet Connection (DSL, Cable)
- Windows 7 or 8 Professional (32 or 64-bit)

MedRx®