



Powerful, small and effective hearing instrument restoration system

ULTRA VAC

The ULTRA VAC system becomes portable with its own custom fitted carrying case.



This durable hard-side case contains a die-cut foam configuration designed especially for the ULTRA VAC and related accessories. Other features include a thick easy grip handle and dual locking system.

Hearing Instrument Restoration System

The Ultra Vac was designed to give dispensers and technicians a functional tool for repairing, restoring and maintaining hearing aids.

MedRx ULTRA VAC

- Pulse / Constant Power Feature
- Durable Vacuum Wand
- Assorted Wand Tips
- Multi-Piece Toolkit
- Pressurized Drying Chamber Cup
- Vacuum Gauge
- Particulate Filter
- Portable with optional Carrying Case

ULTRA VAC

MedRx International

Salzufer 13 / 14
10587 Berlin, Germany
Tel.: +49 30 / 70 71 46-50
Fax: +49 30 / 70 71 46-99
E-Mail: medrx-sales@maico.biz
Web: www.medrx-int.com

Welcome to the New Generation.

Specs

ULTRA VAC

About MedRx

MedRx, Inc. is a U.S. based global manufacturer and developer of advanced computerized diagnostic and hearing instrument fitting technologies, specifically designed for the hearing care professional.

MedRx has created a remarkable new generation of discreet, yet powerful PC-based instrumentation for Audiometry, Real Ear Measurement, Live Speech Mapping, Hearing Instrument Testing & Evaluation and Video Otoscopy.

Standard Accessories

- ULTRA VAC Base Unit
- Vacuum Wand with Hose
- Drying Chamber Cup
- Vacuum Wand Tips
- Cleaning Tool Kit
- ULTRA VAC Case (Optional)

Technical Specifications

POWER REQUIREMENTS

100 VAC, 50-60 Hz
110-120 VAC, 50-60 Hz
220-240 VAC, 50-60 Hz

POWER CONSUMPTION

<100 watts

FUSE

1.25 amp SB 5 x 20 mm
2.50 amp SB 5 x 20 mm

OPERATING TEMPERATURE

0° to 50° C (32° to 120° F)

WEIGHT

Approx. 3.9 kg

DIMENSIONS

(Approx.) 26.7 cm (W)
x 15.9 cm (H)
x 30.5 cm (D)

PARTICULATE FILTER 134S

VACUUM TIPS

#14, #15, 2-#18 and 2-#20 gauge
w/Safety Lock attachment

DRYING CHAMBER

CYCLE TIME

4 minutes with auto shut-off & auto
pressure release

VACUUM/PRESSURE WAND

Continuous use



ULTRA VAC™ is a Trademark of MedRx

Distributed by: